




Gilian® Silica Sampling Kits



The new Gilian Silica Sampling Kits are purpose-built for industrial hygienists needing streamlined, compliant sampling for Respirable Crystalline Silica (RCS).

Available in 3 configurations, each kit includes a Gilian sampling pump (GilAir 3 or GilAir Plus), pre-weighed filter media, and a size selection cyclone – designed to meet NMAM 7500/OSHA silica standards efficiently and accurately.



Kit	Pump & Features	Cyclone/Sampler	Filters	Kit Photo
Basic	GilAir 3	Dorr Oliver 1.7 LPM	Pre-weighed 37mm, 5.0 Micron, 3-pc., 5 filter cassettes	
<i>Calibration jar required, not included. (P/N 7013376)</i>				
Standard	GilAir Plus – includes clock and data logging	Dorr Oliver 1.7 LPM	Pre-weighed 37mm, 5.0 Micron, 3-pc., 5 filter cassettes	
<i>Calibration jar required, not included. (P/N 7013376)</i>				
Advanced	GilAir Plus – includes clock and data logging	GK2.69 Cyclone 4.2 LPM (Connects directly to calibrator, no calibration jar required)	Pre-weighed 37mm, 5.0 Micron, 3-pc., 5 filter cassettes	
<i>Higher flow rate and more reliable results. Captures both thoracic & respirable particles in high and low dust environments.</i>				

Note: Air Flow Calibrator not included with any Kit.

BASIC KIT

Part Number	Description
910-2411-00US-R	Silica Kit, Basic, GilAir3, Dorr-Oliver Cyclone, US Cord
910-2411-00EU-R	Silica Kit, Basic, GilAir3, Dorr-Oliver Cyclone, EU Cord
910-2411-00UK-R	Silica Kit, Basic, GilAir3, Dorr-Oliver Cyclone, UK Cord

STANDARD KIT

Part Number	Description
910-2412-00US-R	Silica Kit, Standard, GilAir Plus, Dorr-Oliver Cyclone, US Cord
910-2412-00EU-R	Silica Kit, Standard, GilAir Plus, Dorr-Oliver Cyclone, EU Cord
910-2412-00UK-R	Silica Kit, Standard, GilAir Plus, Dorr-Oliver Cyclone, UK Cord

ADVANCED KIT

Part Number	Description
910-2413-00US-R	Silica Kit, Advanced, GilAir Plus, GK2.69 Cyclone, US Cord
910-2413-00EU-R	Silica Kit, Advanced, GilAir Plus, GK2.69 Cyclone, EU Cord
910-2413-00UK-R	Silica Kit, Advanced, GilAir Plus, GK2.69 Cyclone, UK Cord

