

NEW!

Gilian® | *MOBILE Pro*



**Operate your Instruments.
Let your Workers Work.**

Remote configuration, operation, and monitoring of devices enabled with Bluetooth®

- » Remotely configure device
- » Remotely operate device
 - Pumps: start, stop, and pause a sample event
 - Calibrator: start a timed validation event, save event record with relevant sample data
- » Remotely monitor pump performance without disturbing the worker or sampling event
- » Remotely review stored / saved events (when device is idle)

GilAir® PLUS



Go-Cal™ Pro



Download on the
App Store



GET IT ON
Google Play



Learn more: [Sensidyne.com/Mobile](https://www.sensidyne.com/Mobile)

SENSIDYNE®
Industrial Health & Safety Instrumentation